



Computer Vision Liver Segmentation with UNet


Installing

```
!pip install nibabel
```

Nibabel 是一個專門用於處理醫學影像的模組
可以用來讀取處理 .nii 或.nii.gz檔案
對於 NIfTI 影像格式支援度非常好。

<https://www.slicer.org/>

可以下載3D Slicer，透過此軟體顯示醫學影像檔案內的影像。



3D Slicer image computing platform

[!\[\]\(1b528d31677da3f213f0576712b16881_img.jpg\) Download](#) [!\[\]\(0f28048ee5f439559f16ca4d39043c34_img.jpg\) Documentation](#)

[!\[\]\(f569861f70431a96430de441727f6f36_img.jpg\) Developers](#) [!\[\]\(8c5d02cd5c755a6d033749bd4db19f3e_img.jpg\) Training](#)

[!\[\]\(b092c9de5ce549057ebaabc0962ed961_img.jpg\) Forum](#) [!\[\]\(53e39e2af75fa6f4f8288923592a70b3_img.jpg\) Twitter](#)

3D Slicer is a **free, open source** and **multi-platform** software package widely used for medical, biomedical, and related imaging research.

Dataset

volume-1.nii	2022/9/20 下午 01:27	NII 檔案	62,977 KB
volume-2.nii	2022/9/20 下午 01:27	NII 檔案	264,705 KB
volume-3.nii	2022/9/20 下午 01:27	NII 檔案	273,409 KB
volume-4.nii	2022/9/20 下午 01:27	NII 檔案	430,593 KB
volume-5.nii	2022/9/20 下午 01:27	NII 檔案	274,945 KB


CT

segmentation-1.nii	2022/9/20 下午 01:25	NII 檔案	31,489 KB
segmentation-2.nii	2022/9/20 下午 01:26	NII 檔案	132,353 KB
segmentation-3.nii	2022/9/20 下午 01:26	NII 檔案	136,705 KB
segmentation-4.nii	2022/9/20 下午 01:26	NII 檔案	430,593 KB
segmentation-5.nii	2022/9/20 下午 01:26	NII 檔案	137,473 KB


Mask

Dataset

CT



(a) coronal view






(b) axial view



(c) sagittal view

Mask




Data


1. 讀取Npy檔，上述醫學檔案已處理好儲存在一起。

```
total_patch = np.load("./drive/MyDrive/Unet/trainNpy/total_patch.npy")
total_mask = np.load("./drive/MyDrive/Unet/trainNpy/total_mask.npy")

total_patch.shape
(15318, 64, 64, 1)
```



Unet Model



LOSS Function

Weighted binary crossentropy 用於二分類

$$\text{BCELoss}(O, T) = -\frac{1}{n} \sum_i ((T[i] * \log(O[i])) + (1 - T[i]) * \log(1 - O[i]))$$

True Prob of True Prob of
label Positive label Negative
 Class Class

```
#loss function
def weighted_binary_crossentropy(y_true, y_pred):
    y_pred = tf.clip_by_value(y_pred, 10e-8, 1.-10e-8) #把值壓縮到(min,max)之間，小於min的值變成min，大於max的值變成max
    loss = - (y_true * K.log(y_pred) * 0.90 + (1 - y_true) * K.log(1 - y_pred) * 0.10)


    return K.mean(loss)
```

因為正樣本(擁有肝臟的Mask)較少，因此要給予較大的權重。

LOSS Function

Dice Coefficient Loss 計算樣本相似度(數值越大即重疊率越高)

```
smooth = 1. #用於防止分母為0
def dice_coef(y_true, y_pred):
    y_true_f = K.flatten(y_true) #攤平成一維
    y_pred_f = K.flatten(y_pred)
    intersection = K.sum(y_true_f * y_pred_f)
    return (2. * intersection + smooth) / (K.sum(y_true_f) + K.sum(y_pred_f) + smooth)
```



Training

1. 設定訓練過程需要的優化器、損失函數和評估標準

- optimizer : 優化器，梯度下降，使Loss可以越小越好
- loss : 實際值和預測值的差距，是優化器要指引的對象
- metrics : 與Loss相似，但結果並不會回饋到訓練過程中

```
adam = Adam(lr = 0.0001)
model.compile(optimizer = adam, loss = weighted_binary_crossentropy, metrics = [dice_coef])
```

2. 設定參數，開始訓練

- total_patch : 訓練資料
- total_mask : 訓練資料的答案
- batch_size : 一次訓練的樣本數
- epoch : 設定訓練整個資料集的次數

```
model.fit(total_patch, total_mask, batch_size = 128, epochs = 30)
```

Training




```
Epoch 41/50
120/120 [=====] - 27s 223ms/step - loss: 0.0245 - dice_coef: 0.9090
Epoch 42/50
120/120 [=====] - 27s 223ms/step - loss: 0.0238 - dice_coef: 0.9099
Epoch 43/50
120/120 [=====] - 27s 223ms/step - loss: 0.0266 - dice_coef: 0.9043
Epoch 44/50
120/120 [=====] - 27s 223ms/step - loss: 0.0298 - dice_coef: 0.8990
Epoch 45/50
120/120 [=====] - 27s 223ms/step - loss: 0.0247 - dice_coef: 0.9097
Epoch 46/50
120/120 [=====] - 27s 223ms/step - loss: 0.0346 - dice_coef: 0.8993
Epoch 47/50
120/120 [=====] - 27s 222ms/step - loss: 0.0263 - dice_coef: 0.9010
Epoch 48/50
120/120 [=====] - 27s 223ms/step - loss: 0.0228 - dice_coef: 0.9118
Epoch 49/50
120/120 [=====] - 27s 222ms/step - loss: 0.0217 - dice_coef: 0.9170
Epoch 50/50
120/120 [=====] - 27s 222ms/step - loss: 0.0201 - dice_coef: 0.9177
<keras.callbacks.History at 0x7f3e1c094a60>
```

Saving Model

```
#儲存模型
model_json = model.to_json()
with open("./model_json_final.json", "w") as json_file:
    json_file.write(model_json)


# serialize weights to HDF5
model.save_weights("./model_weights_final.h5")
print("Saved model to disk")
```



```
# load weights into new model
json_file = open("./model_json_final.json", 'r')
loaded_model_json = json_file.read()
json_file.close()

loaded_model = model_from_json(loaded_model_json)
loaded_model.load_weights("./model_weights_final.h5")
print("Loaded model from disk")
```

Prediction






```
8/8 [=====] - 0s 12ms/step
```

```
Slice to Patch Shape: (225, 64, 64, 1)
```

```
Prediction Shape: (512, 512, 1)
```

```
Dice coef: tf.Tensor(0.8556331152982525, shape=(), dtype=float64)
```

Prediction



8/8 [=====] - 0s 12ms/step

Slice to Patch Shape: (225, 64, 64, 1)

Prediction Shape: (512, 512, 1)

Dice coef: tf.Tensor(0.8408137781547451, shape=(), dtype=float64)



Thank you